

## Technical Data Sheet for Glass Microfibre Syringe Filters (25mm Diameter)

These are high-viscosity filters with a pre-filter to enable filtration of aqueous based biorelevant dissolution samples.

## **Technical Information**

25mm diameter, 0.45 micrometer pore size, GMF with glass pre-filter
GMF25
Filtration of dissolution samples where sampling volume is >5mL
Dissolution with both fasted and fed state biorelevant media
Yes
Polypropylene
30mm x 21mm
>5mL
3g
1.4mL with air purge 250µL
72
Not sterile but can be autoclaved at 121°C 11 bar for 20 minutes
Inlet - Female Luer Lock (FLL) Outlet - Male Slip Luer (MSL)
Stacked within plastic bag
Bead retention and burst pressure
84212900
25°C 15°C

For technical help, contact <a href="mailto:support@biorelevant.com">support@biorelevant.com</a>